

# Syringe Filter Chemical Compatibility

Chemical compatibility is a critical consideration when selecting the proper sample prep syringe filter for your application. This chart outlines the chemical compatibility of the most common syringe filters.

| Syringe Filter Chemical Compatibility | MEMBRANE                   |                                |                        |   |                        |            |                  | HOUSING |      |
|---------------------------------------|----------------------------|--------------------------------|------------------------|---|------------------------|------------|------------------|---------|------|
|                                       | Regenerated Cellulose (RC) | Polytetrafluoroethylene (PTFE) | Cellulose Acetate (CA) | Cellulose Acetate + Glass Fiber (CA + GF) | Polyethersulfone (PES) | Nylon (NY) | Glass Fiber (GF) | MBS     | PP   |
| Filter Housing                        | PP                         | PP                             | MBS                    | MBS                                       | MBS                    | PP         | MBS              | MBS     | PP   |
| <b>Sterilization</b>                  |                            |                                |                        |   |                        |            |                  |         |      |
| Ethylene oxide                        | ++                         | ++                             | ++                     | ++  | ++                     | ++         | ++               | ++      | ++   |
| Gamma irradiation                     | -                          | -                              | ++                     | ++  | ++                     | -          | ++               | ++      | -    |
| Autoclaving 121 °C, 30 min            | ++                         | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| <b>Solvents</b>                       |                            |                                |                        |   |                        |            |                  |         |      |
| Acetone                               | ++                         | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| Acetonitrile                          | ++                         | ++                             | -                      | -   | -                      | n.a.       | -                | -       | ++   |
| Benzene                               | ++                         | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| Benzyl alcohol                        | +                          | +                              | -                      | -   | -                      | +          | -                | -       | +    |
| n-Butyl acetate                       | ++                         | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| n-Butanol                             | ++                         | ++                             | +                      | +   | +                      | ++         | ++               | ++      | ++   |
| Carbon tetrachloride                  | -                          | -                              | -                      | -   | -                      | -          | -                | -       | -    |
| Cellosolve                            | -                          | -                              | -                      | -   | -                      | -          | -                | -       | -    |
| Chloroform                            | ++                         | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| Cyclohexane                           | +                          | +                              | +                      | +   | -                      | +          | +                | +       | +    |
| Cyclohexanone                         | +                          | +                              | -                      | -   | -                      | +          | -                | -       | +    |
| Diethylacetamide                      | ++                         | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| Diethyl ether                         | ++                         | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| Dimethyl formamide                    | +                          | +                              | -                      | -   | -                      | +          | -                | -       | +    |
| Dimethylsulfoxide                     | ++                         | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| Dioxane                               | ++                         | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| Ethanol, 98%                          | +                          | +                              | -                      | -   | -                      | +          | -                | -       | +    |
| Ethyl acetate                         | +                          | +                              | -                      | -   | -                      | +          | -                | -       | +    |
| Ethylene glycol                       | ++                         | ++                             | +                      | +   | ++                     | ++         | ++               | ++      | ++   |
| Formamide                             | +                          | ++                             | -                      | -   | ++                     | ++         | ++               | ++      | ++   |
| Glycerin                              | +                          | +                              | +                      | +   | +                      | +          | +                | +       | +    |
| n-Heptane                             | ++                         | ++                             | +                      | +   | +                      | ++         | +                | +       | ++   |
| n-Hexane                              | +                          | +                              | +                      | +   | +                      | +          | +                | +       | +    |
| Isobutanol                            | -                          | -                              | +                      | +   | ++                     | -          | ++               | ++      | -    |
| Isopropanol                           | ++                         | ++                             | -                      | +   | -                      | ++         | -                | -       | ++   |
| Isopropyl acetate                     | ++                         | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| Methanol, 98%                         | +                          | +                              | -                      | +   | +                      | +          | ++               | ++      | +    |
| Methyl acetate                        | +                          | +                              | -                      | -   | -                      | +          | -                | -       | +    |
| Methylene chloride                    | ++                         | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| Methyl ethyl ketone                   | +                          | +                              | -                      | -   | -                      | +          | -                | -       | +    |
| Methyl isobutyl ketone                | +                          | +                              | -                      | -   | -                      | +          | -                | -       | +    |
| Monochlorobenzene                     | +                          | +                              | -                      | -   | -                      | +          | -                | -       | +    |
| Nitrobenzene                          | +                          | +                              | -                      | +   | -                      | +          | -                | -       | +    |
| n-Pentane                             | ++                         | ++                             | +                      | +   | +                      | ++         | +                | +       | ++   |
| Perchloroethylene                     | ++                         | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| Pyridine                              | ++                         | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| Tetrahydrofuran                       | ++                         | ++                             | -                      | ++  | -                      | ++         | -                | -       | ++   |
| Toluene                               | ++                         | ++                             | -                      | ++  | -                      | ++         | -                | -       | ++   |
| Trichloroethane                       | n.a.                       | n.a.                           | -                      | -   | -                      | n.a.       | -                | -       | n.a. |
| Trichloroethylene                     | ++                         | ++                             | -                      | ++  | -                      | ++         | -                | -       | ++   |
| Xylene                                | +                          | +                              | -                      | -   | -                      | +          | -                | -       | +    |
| <b>Acids</b>                          |                            |                                |                        |   |                        |            |                  |         |      |
| Acetic acid, 25%                      | +                          | +                              | -                      | -   | -                      | -          | -                | -       | +    |
| Acetic acid, 80%                      | +                          | +                              | -                      | -   | -                      | -          | -                | -       | +    |
| Hydrofluoric acid, 25%                | +                          | +                              | -                      | -   | +                      | -          | +                | +       | +    |
| Hydrofluoric acid, 50%                | +                          | +                              | -                      | -   | +                      | -          | +                | +       | +    |
| Hydrochloric acid, 15%                | -                          | +                              | +                      | +   | +                      | -          | +                | +       | +    |
| Hydrochloric acid, 20%                | -                          | +                              | -                      | -   | +                      | -          | +                | +       | +    |
| Nitric acid, 30%                      | -                          | +                              | -                      | -   | +                      | -          | +                | +       | +    |
| Nitric acid, conc.                    | -                          | -                              | -                      | -   | -                      | -          | -                | -       | -    |
| Perchloric acid, 25%                  | -                          | +                              | -                      | -   | -                      | -          | n.a.             | n.a.    | +    |
| Phosphoric acid, 1%                   | -                          | +                              | +                      | +   | +                      | -          | +                | +       | +    |
| Phosphoric acid, 86%                  | -                          | +                              | +                      | +   | +                      | -          | +                | +       | +    |
| Sulfuric acid, 25%                    | +                          | ++                             | -                      | -   | +                      | -          | +                | +       | ++   |
| Sulfuric acid, 98%                    | -                          | +                              | -                      | -   | -                      | -          | -                | -       | +    |
| Trichloroacetic acid, 25%             | +                          | +                              | -                      | -   | -                      | -          | -                | -       | +    |
| <b>Bases</b>                          |                            |                                |                        |   |                        |            |                  |         |      |
| Ammonia, 1N                           | +                          | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| Ammonium hydroxide, 25%               | +                          | +                              | -                      | -   | -                      | +          | -                | -       | +    |
| Potassium hydroxide, 32%              | -                          | ++                             | -                      | -   | -                      | +          | -                | -       | ++   |
| Sodium hydroxide, 32%                 | -                          | +                              | -                      | -   | -                      | +          | -                | -       | +    |
| Sodium hydroxide, 1N                  | +                          | ++                             | -                      | -   | -                      | ++         | -                | -       | ++   |
| <b>Aqueous Solutions</b>              |                            |                                |                        |   |                        |            |                  |         |      |
| Formalin, 30%                         | +                          | +                              | +                      | +   | +                      | +          | +                | +       | +    |
| Sodium hypochlorite, 5%               | -                          | +                              | -                      | -   | +                      | -          | +                | +       | +    |
| Hydrogen peroxide, 35%                | -                          | ++                             | -                      | -   | +                      | -          | +                | +       | ++   |
| <b>pH Range</b>                       | pH 1-14                    |                                |                        |   |                        |            |                  | -       | ++   |
|                                       | pH 1-13                    |                                |                        |   |                        |            |                  | -       | ++   |
|                                       | pH 3-14                    |                                |                        |   |                        |            |                  | -       | ++   |
|                                       | pH 3-12                    |                                |                        |   |                        |            |                  | +       | ++   |
|                                       | pH 4-8                     |                                |                        |   |                        |            |                  | ++      | ++   |

#### Legend

Compatible: ++  
 Limited compatibility: +  
 Not compatible: -  
 MBS: Methacrylate Butadiene Styrene  
 PP: Polypropylene  
 n.a.: Not analyzed

Contact time: 24 hours at 20 °C  
 Chemical compatibilities can be influenced by various factors. Therefore, we recommend that you confirm compatibility with the liquid you want to filter by performing a trial filtration run before you start your actual filtration. Both membrane & housing compatibility need to be considered together.

6008.L



Trademarks

Phenex is a trademark of Phenomenex, Inc. Teflon is a registered trademark of E.I. du Pont de Nemours and Co.

Disclaimer

Subject to Phenomenex standard Terms and Conditions which may be viewed at [www.Phenomenex.com/TermsAndConditions](http://www.Phenomenex.com/TermsAndConditions)

© 2008 Phenomenex, Inc. All rights reserved.

Phenomenex | Tel: 310-212-0555 | Fax: 310-328-7768 | Email: [info@phenomenex.com](mailto:info@phenomenex.com) | Web: [www.phenomenex.com](http://www.phenomenex.com)

