


INSTALLATION INSTRUCTIONS

The SecurityLINK UHPLC fingertight fitting system simplifies your system connection process and provides consistent performance with torque limiting technology that prevents column damaging overtightening. No tools required for installation again!

- 1 Insert SecurityLINK UHPLC fitting into column port



- 2 Finger tighten until first "CLICK" feedback is received.* Connection Complete!



Tip

The "CLICK" Feedback indicates the SecurityLINK Connection is Secure!
This Prevents Overtightening & Saves Your Column.

ORDERING INFORMATION

PEEKsil™ Double Sided Fittings



ID (µm)	LENGTH (mm)	Fitting Size Top (in)	Fitting Size Bottom (in)	Part No
25	100	1/16	1/16	AJ1-2111
25	150	1/16	1/16	AJ1-2121
25	250	1/16	1/16	AJ1-2141
25	300	1/16	1/16	AJ1-2151
25	500	1/16	1/16	AJ1-2171
25	750	1/16	1/16	AJ1-2191
25	1000	1/16	1/16	AJ1-21A1
50	100	1/16	1/16	AJ1-2211
50	150	1/16	1/16	AJ1-2221
50	200	1/16	1/16	AJ1-2231
50	250	1/16	1/16	AJ1-2241
50	300	1/16	1/16	AJ1-2251
50	500	1/16	1/16	AJ1-2271
50	750	1/16	1/16	AJ1-2291
50	1000	1/16	1/16	AJ1-22A1
75	150	1/16	1/16	AJ1-2321
75	250	1/16	1/16	AJ1-2341
75	500	1/16	1/16	AJ1-2371
75	1000	1/16	1/16	AJ1-23A1
100	100	1/16	1/16	AJ1-2411
100	150	1/16	1/16	AJ1-2421
100	250	1/16	1/16	AJ1-2441
100	500	1/16	1/16	AJ1-2471
100	1000	1/16	1/16	AJ1-24A1

Stainless Steel Double Sided Fittings




ID (µm)	LENGTH (mm)	Fitting Size Top (in)	Fitting Size Bottom (in)	Part No
100	150	1/16	1/16	AJ1-1421
100	250	1/16	1/16	AJ1-1441
100	350	1/16	1/16	AJ1-1461
100	500	1/16	1/16	AJ1-1471
100	600	1/16	1/16	AJ1-1481
125	150	1/16	1/16	AJ1-1521
125	250	1/16	1/16	AJ1-1541
125	350	1/16	1/16	AJ1-1561
125	500	1/16	1/16	AJ1-1571
125	600	1/16	1/16	AJ1-1581
254	150	1/16	1/16	AJ1-1621
254	250	1/16	1/16	AJ1-1641
254	350	1/16	1/16	AJ1-1661
254	500	1/16	1/16	AJ1-1671
254	600	1/16	1/16	AJ1-1681

PEEK-Lined Stainless Steel Double Sided Fittings



ID (µm)	LENGTH (mm)	Fitting Size Top (in)	Fitting Size Bottom (in)	Part No
25	150	1/16	1/16	AJ1-3121
25	250	1/16	1/16	AJ1-3141
25	350	1/16	1/16	AJ1-3161
25	500	1/16	1/16	AJ1-3171
25	600	1/16	1/16	AJ1-3181
50	150	1/16	1/16	AJ1-3221
50	250	1/16	1/16	AJ1-3241
50	350	1/16	1/16	AJ1-3261
50	500	1/16	1/16	AJ1-3271
50	600	1/16	1/16	AJ1-3281
75	150	1/16	1/16	AJ1-3321
75	250	1/16	1/16	AJ1-3341
75	350	1/16	1/16	AJ1-3361
75	500	1/16	1/16	AJ1-3371
75	600	1/16	1/16	AJ1-3381
100	150	1/16	1/16	AJ1-3421
100	250	1/16	1/16	AJ1-3441
100	350	1/16	1/16	AJ1-3461
100	500	1/16	1/16	AJ1-3471
100	600	1/16	1/16	AJ1-3481
75	250	1/16	1/16	AJ1-2341
75	500	1/16	1/16	AJ1-2371
75	1000	1/16	1/16	AJ1-23A1
100	100	1/16	1/16	AJ1-2411
100	150	1/16	1/16	AJ1-2421
100	250	1/16	1/16	AJ1-2441
100	500	1/16	1/16	AJ1-2471
100	1000	1/16	1/16	AJ1-24A1

PEEKsil Single Sided Fittings



ID (µm)	LENGTH (mm)	Fitting Size Top (in)	Fitting Size Bottom (in)	Part No
50	150	1/16	None	AJ1-2224
50	500	1/16	None	AJ1-2274
50	750	1/16	None	AJ1-2294
50	1000	1/16	None	AJ1-22A4

Reorder Tip:

SecurityLINK tubing material includes a sleeve that provides: ID, length, and part number information.



Terms and Conditions

Subject to Phenomenex standard terms and conditions which can be viewed at www.phenomenex.com/TermsandConditions

Trademarks

SecurityLINK is a trademark of Phenomenex. PEEKsil is a trademark of SGE International Pty. Ltd.

©2018 Phenomenex, Inc. All rights reserved.



If SecurityLINK does not perform as well or better than your current UHPLC fitting system of similar type, dimensions, and material, return the product with comparative data within 45 days for a FULL REFUND.